

thermage®
FLX
RF system

GIVE YOUR COLLAGEN A WAKE-UP CALL

**OVER
3.3 MILLION**
TREATMENTS
PERFORMED

MODEL AND ACTUAL PATIENT

GET BACK TO THE REAL YOU

Thermage FLX is a non-invasive treatment that uses monopolar radiofrequency technology to smooth, tighten, and contour the eyes, face, and body for an overall younger-looking appearance. Thermage FLX can be used to treat both genders of all ages and all skin types, delivering lasting results in a single treatment.



EASY TO NAVIGATE

Large touchscreen with user friendly interface.



VERSATILE & USER-FRIENDLY

One handpiece to treat the face, eyes and body. Even easier to use with a foot pedal.



PATIENT COMFORT

Vibration available for both face and body treatments.



PATENTED ACCUREP TECHNOLOGY

Optimised treatment algorithm with site specific real-time tuning. Faster treatment time and a single handpiece for all treatment areas.



MULTIPLE TREATMENT AREAS

Designed to improve wrinkles and loose skin on the face, eyelids and body.



ENHANCED USER EXPERIENCE

Lighter, ergonomic handpiece with multi-directional vibration.



Patients reported being clearly satisfied with the Thermage FLX Treatment¹

9.4

AVERAGE SCORE




1. Data on file, Social Medical


SMOOTHER, TIGHTER, MORE CONTOURED SKIN




THERMAGE FLX HELPS REDUCE THE VISIBLE SIGNS OF AGING IN MINUTES




As we age, the collagen that keeps skin firm breaks down resulting in sagging skin, fine lines, and wrinkles.



The latest generation of Thermage® uses radiofrequency technology to heat the deep, collagen-rich layers of the skin, while the tip vibrates and cools the surface to help aid in patient comfort.




The applied heat separates the water molecules from the collagen fibrils, causing an immediate contraction, resulting in skin tightening.



Over the next few months, the secondary healing response continues as collagen is deposited and remodeled.



This process tightens the skin further over time.



You will notice smoother, tighter skin just after one treatment with continued improvement for up to 6 months.

CUSTOMIZE YOUR TREATMENTS WITH THE RIGHT TIPS

Enhance your treatment outcomes by using the appropriate Thermage FLX tips for each treatment area. Uses precise heating to address wrinkles. Total tip and Body tip offers vibration to help aid in patient comfort.



Total Tip
4.0 cm²

for treatments on- and off-face



Eye Tip
0.25 cm²

for treatments on the eyelids only



Total Tip*
3.0 cm²

For peri-orbital area



Body Tip
16 cm²

for treatments on body areas only

RESULTS THAT MATTER



Before



1 Month Post Treatment



Before



2 Months Post Treatment



Before



3 Months Post Treatment

Before and After photos have not been retouched and represent patients who have had significant results. Individual results may vary based on the aggressiveness of treatment and skin condition.

SOLTA SUPPORT TO DRIVE SUCCESS



Training Support

Solta Academy is an online educational platform that provides opportunities to physicians and their staff to build on their clinical knowledge and gain insights at their own pace and time to help deliver exceptional patient care.



Practice Support

Offering assistance with business planning, marketing strategies and support. Includes placement on Solta Medical's clinic finder site, customised bespoke marketing materials alongside a full suite of ongoing complimentary digital and print resources, ensuring your marketing efforts are always effective.



Product Support

Standard 1-year manufacturer warranties with access to responsive Product Support Team that includes field service engineers for assistance in minimising business disruption.



SOLTAMEDICAL

This product is not available for purchase by the general public. Except as otherwise indicated, all product names, slogans, and other marks are trademarks of the Bausch Health group of companies. © 2025 Solta Medical, a division of Bausch Health Companies Inc. All rights reserved. Distributed in Australia by Solta Medical Australia Pty Ltd ABN 71 652 701 767. Suite 301, Level 3, 12 Help St, Chatswood, NSW 2067. (02) 8311 0305. Distributed in New Zealand by Radiant Health Ltd NZBN 9429032273090. Pier 21 Centre, Level 3, 11 Westhaven Drive, Auckland Central 1010. 0508 35394. THR.0018.AU.25 G247- 106600



FOLLOW US ON

Solta Medical Australia

@Thermageau

www.thermage.com.au

Thermage AU